

Comparison of a highly sensitive, tumor-informed Illumina whole genome sequencing research assay and a tumor-informed bespoke whole exome sequencing assay for molecular residual disease detection in CheckMate 77T (NCT 04025879)



8591



Jonathan Beer<sup>1</sup>, Khai Luong<sup>2</sup>, Jeff Tsai<sup>2</sup>, Jacob Gibson<sup>2</sup>, Jonathan Baden<sup>1</sup>, Sven Bilke<sup>2</sup>, Kate Broadbent<sup>2</sup>, Yuan Ding<sup>2</sup>, Christine Eng<sup>2</sup>, Yu-Han Hung<sup>1</sup>, Artem Movsesyan<sup>2</sup>, Jaclyn Neely<sup>1</sup>, Mark Stern<sup>1</sup>, Weixin Wu<sup>2</sup>, David Yu<sup>1</sup>, Yunjiao Zhu<sup>2</sup>, Traci Pawlowski<sup>2</sup>, Li Liu<sup>2\*</sup>

<sup>1</sup>Bristol Myers Squibb; <sup>2</sup>Illumina; \*corresponding author

Background

- Brief overview of CheckMate 77T and clinical context:
  - Phase III trial of nivolumab (NIVO) vs. placebo (PBO) as perioperative therapy
  - 461 patients with resectable stage IIA-IIIB NSCLC
  - Clinically meaningful improvements in event-free survival (EFS) and pathologic complete response (pCR) with NIVO observed<sup>1</sup>
- Prior ctDNA findings<sup>2-4</sup>:
  - ctDNA clearance during neoadjuvant period: 66% (NIVO) vs. 38% (PBO)
  - Reported from 130 patients using a tumor-informed bespoke panel (TIBP) WES MRD assay
- Gap / Rationale:
  - TIBP WES MRD assays have high operational complexity and long turnaround time (TAT)
  - Hypothesis: A WGS-based MRD assay may offer improved sensitivity and simplified workflow
- Study Objective: Compare ctDNA clearance rates and performance metrics between the Illumina WGS MRD assay and the TIBP WES MRD assay in a subset of CheckMate 77T participants

Methods

- Study Design / Patient Population:
  - Substudy of CheckMate 77T
  - 13% (60/461) of total study population, all with TIBP WES MRD data previously reported
  - WGS MRD from residual tumor FPPE tissue or isolated DNA, PBMC germline normal, cfDNA from plasma
- Illumina WGS MRD Assay Workflow:
  - Subtracting WGS germline normal from WGS of tumor → tumor-specific genomic fingerprint
  - Bioinformatic monitoring of plasma timepoints via low-pass WGS (lpWGS)
  - Library prep: Illumina WGS Oncology Prep (research use only)
  - Analysis: DRAGEN WGS MRD pipeline
  - Analytical sensitivity and specificity: LoD95 as low as 10 ppm at 99.5% specificity from 5 ng
- Tumor Informed Bespoke Panel Workflow:
  - TIBP WES MRD assay (used in prior Provencio et al. analysis)
- Key Analyses:
  - ctDNA detection rates at NEO C1D1 (prior to neoadjuvant treatment)
  - ctDNA clearance rates during the neoadjuvant period (NEO C1D1 → NEO END) by treatment arm (NIVO vs. PBO)
  - Concordance metrics: PPA, NPA, OPA (TIBP WES as reference standard)
  - Head-to-head comparison of operational complexity and TAT

Results

Table 1. Comparison of MRD assay characteristics

| Parameter                     | Illumina WGS       | TIBP WES           |
|-------------------------------|--------------------|--------------------|
| Patient-specific panel design | Not required       | Required           |
| cfDNA input requirement       | 2-5 ng             | Typically 10-60 ng |
| Workflow                      | Simple streamlined | Complex multi-step |
| TAT*                          | 1-2 weeks          | 4-6 weeks          |
| Distributable kit             | Yes                | No                 |

\*TAT (turnaround time) for tumor tissue, PBMC germline, and first plasma timepoint results

Table 2. Samples included in the analysis

| ARM             | Timepoint | Patient count |
|-----------------|-----------|---------------|
| NIVO+CHEMO_NIVO | NEO C1D1  | 27            |
| NIVO+CHEMO_NIVO | NEO END   | 26            |
| PBO+CHEMO_PBO   | NEO C1D1  | 30            |
| PBO+CHEMO_PBO   | NEO END   | 28            |

- TAT is 3-4x shorter and cfDNA input requirement is 5-12x lower for Illumina WGS (Table 1, Figure 1)
- Only samples that meet the following two criteria were included in this analysis (Table 2).
  - Samples collected during neoadjuvant phase (NEO C1D1, NEO END)
  - Samples that have both TIBP WES and Illumina WGS readout
- Filtering yields 111 samples from 60 unique patients

Figure 1. Schematic diagram of MRD workflows




Table 3. ctDNA Detection at NEO C1D1 and NEO END

| Assay                           | NEO C1D1 (n=57) | NEO END (n=42) | Overall (n=111) |
|---------------------------------|-----------------|----------------|-----------------|
| Illumina WGS MRD (this study)   | 91%             | 48%            | 70%             |
| TIBP WES MRD (Provencio et al.) | 86%             | 41%            | 64%             |

Table 4. Confusion Matrix for TIBP WES (rows) vs Illumina WGS (columns)

|         | ABSENT | PRESENT |    |
|---------|--------|---------|----|
| ABSENT  | 32     | 8       | 40 |
| PRESENT | 1      | 70      | 71 |
|         | 33     | 78      |    |


Table 5. Concordance with TIBP WES (reference) at NEO C1D1 and NEO END

|     | NEO C1D1 | NEO END | Overall |
|-----|----------|---------|---------|
| PPA | 100%     | 95%     | 99%     |
| NPA | 63%      | 84%     | 80%     |
| OPA | 95%      | 89%     | 92%     |

PPA: positive percent agreement; NPA: negative percent agreement; OPA: overall percent agreement

- ctDNA detection rates prior to neoadjuvant treatment (NEO C1D1) were greater in this study than previously reported, 91% for Illumina WGS vs 86% for TIBP WES (Table 3)
- 8 out of 9 disagreement cases (3 at NEO C1D1 and 6 at NEO END) were samples classified as “ctDNA present” by Illumina WGS but “ctDNA absent” by TIBP WES (Table 4)
  - Samples collected during neoadjuvant phase (NEO C1D1, NEO END)
  - Samples that have both TIBP WES and Illumina WGS readout
- A high overall concordance rate (92% OPA) was observed between the TIBP WES and Illumina WGS platforms. The concordance rate prior to neoadjuvant treatment was 95% and 89% at the end of neoadjuvant treatment (Table 5)

Figure 3. EFS by ctDNA clearance during neoadjuvant phase



- Both platforms showed that ctDNA clearance is associated with better EFS. Hazard Ratio between clearance and non-clearance outcomes is 0.25 with Illumina WGS MRD vs 0.33 with TIBP WES (Figure 3)

Conclusions

- The Illumina WGS MRD assay demonstrated high concordance with tumor-informed WES MRD for MRD-positive detection (PPA = 99%)
- Illumina WGS MRD demonstrated improved analytical sensitivity at baseline (91% vs. 86% detection)
- Lower Hazard Ratio between clearance and non-clearance outcomes with Illumina WGS MRD (0.25 vs. 0.33)
- Significantly reduced TAT with Illumina WGS MRD (1-2 weeks vs. 4-6 weeks)
- Lower cfDNA input requirements with Illumina WGS MRD (2-5 ng vs. 10-60 ng)
- Decreased operational complexity of Illumina WGS MRD — eliminates need for patient-specific panel design, manufacturing, and validation
- Illumina’s WGS-based MRD may represent a distributable, clinically scalable and operationally practical approach for ctDNA monitoring in resectable NSCLC

References

- Cascone et al. Perioperative Nivolumab in Resectable Lung Cancer. NEJM, 2024
- Provencio et al. LBA50, ESMO 2024
- Cascone et al. LBA8010, ASCO 2025
- Cascone et al. Presentation number CT015, AACR 2026

Acknowledgments & Disclosures

- The patients and families who made this study possible
- The clinical study teams who participated in CheckMate 77T
- Bristol Myers Squibb (Princeton, NJ, USA) and Ono Pharmaceutical Company Ltd. (Osaka, Japan)
- The study was supported by Bristol Myers Squibb; WGS MRD analysis was supported by Illumina
- Illumina WGS Oncology Prep is for research use only

Scientific Content on Demand



QR codes are valid for 90 days after the congress presentation date.